

US Distributor of Magstim Magnetic Stimulators



JALI MEDICAL

introduces

magstim®

Air Film Coil



NOW YOU CAN DO MORE

The Magstim® Air Film™ Coil is the first in a new generation of stimulating coils which allow users to stimulate for extended periods of time. Using managed ambient air flow and integrated temperature-regulated fan technology, the Air Film Coil is cooled without the need of an external cooling supply or water reservoir. The novel air cushion used prevents the coil from ever contacting the patient, resulting in uncompromised patient safety.

INTENDED USE: Stimulation of Peripheral Nerves.

Air Coil Film specifications

Stimulating Versatility

The Air Film Coil has been created to work in combination with Magstim Rapid² stimulators, allowing operators to utilize the full capabilities of the range.

Totally Flexible

This unique Air Film Coil is an integral part of a package which also includes a highly flexible and ergonomically designed positioning arm, which facilitates extended periods of stimulation. Fully articulated, with friction joints, the Air Film Coil arm makes handling and positioning simple.

A placebo version is also available for investigational use. The Air Film Placebo Coil looks, sounds and feels the same as a live coil both to the subject and operator, but does not deliver deep stimulation.

BENEFITS:

- Delivers up to of 3,000 (Biphasic waveform) over a 30 minute period and remains within allowed temperature limits
- Extends the useful range of stimulating protocols available to researchers
- Indefinite running time with certain protocols
- Uncompromised patient safety
- Quiet running
- Easy to use, ergonomic design

INTENDED USE:
Stimulation of Peripheral Nerves.

JALI MEDICAL, INC. Inquiries: 888-959-JALI

✉ 500 West Cummings Park, Suite 4950
Woburn, MA 01801

☎ 781-376-0808

✉ inquiries@jalimedical.com

🌐 www.jalimedical.com

October 2008

