

US Distributor of Magstim Magnetic Stimulators



# JALI MEDICAL

*introduces*

**magstim**<sup>®</sup>

## *BiStim<sup>2</sup>* *Magnetic Stimulator*



PROGRAMMABLE INTERVAL PAIRED PULSE STIMULATION

WITH THE ADDED FLEXIBILITY OF A COMBINED PULSE

BiStim<sup>2</sup> provides an effective upgrade path that allows users to take full advantage of the Magstim 200<sup>2</sup> modularity. The potential to combine two units offers fully programmable paired pulse stimulation through a single stimulating coil. The BiStim<sup>2</sup> is capable of producing close interval pulse pairs.



CAUTION: INVESTIGATIONAL DEVICE. Federal (or US) law limits device to investigational use. For information regarding research, please visit [www.jalimedical.com](http://www.jalimedical.com).

# Bistim<sup>2</sup> specifications

BiStim<sup>2</sup> has the added advantage of the ability to summate the two single pulses provided by Magstim stimulators to produce a single High Power pulse.

The interval between stimuli is adjustable using either the integral stimulator controls, remote control coil or externally controlled via triggering software, offering complete user flexibility.

## Two inter pulse spacing options offering maximum controllability:

- 1ms - 999ms in 1ms increments
- 1.0ms - 99.9ms in 0.1ms increments

BiStim<sup>2</sup> remote control coil options allow complete 'single handed' system control including:

- amplitude manipulation of both stimulators
- inter pulse space selection.

## Output Type:

Monophasic Twin Pulse

## Maximum Power as % of Magstim 200<sup>2</sup>:

113% \*

## Maximum Frequency at Maximum Power:

0.2Hz

## Upper Frequency:

0.2Hz

## Maximum Burst Frequency:

1000Hz

## Minimum Pulse Interval:

1 msec

\* in Single High Pulse mode



**CAUTION: INVESTIGATIONAL DEVICE.** Federal (or US) law limits device to investigational use.  
For information regarding research, please visit [www.jalimedical.com](http://www.jalimedical.com).

JALI MEDICAL, INC. Inquiries: 888-959-JALI

✉ 500 West Cummings Park, Suite 4950

Woburn, MA 01801

☎ 781-376-0808

✉ [inquiries@jalimedical.com](mailto:inquiries@jalimedical.com)

🌐 [www.jalimedical.com](http://www.jalimedical.com)

October 2008



US Distributor of Magstim Magnetic Stimulators

© 2008 Jali Medical, Inc.