

US Distributor of Magstim Magnetic Stimulators



JALI MEDICAL

introduces

magstim[®]



Magstim Lab Magnetic Stimulator

CREATING THE FUTURE THROUGH INNOVATION

Magstim Lab biphasic magnetic stimulator draws upon the years of knowledge and expertise gathered in the research environment. Its innovative and compact design allows for the delivery of repetitive stimulation at a fraction of the cost of other repetitive stimulators. The device is simple to operate and has been designed specifically for use with small coils in the laboratory setting.



CAUTION: Device for investigational use in laboratory animals or other test that does not involve human subjects.

Magstim Lab specifications

Power Input:

110 - 120V, 50 - 60Hz

Pulse Characteristics:

Biphasic

Repetition Rate:

0 to 30% - 10 Hz

31 to 80% - 5 Hz

81 to 100% - 3 Hz

Capacitor Life Expectancy:

2×10^6 discharges at 70% power level

8×10^5 discharges at 80% power level

4×10^5 discharges at 90% power level

2×10^5 discharges at 100% power level

UI Controller:

Easy to read digital display of the power level, coil temperature and system status

Approx. Dimensions:

18.1" Wide x 6.3" High x 14.8" Deep

Approx. Weight:

45 LB



JALI MEDICAL, INC. Inquiries: 888-959-JALI

✉ 500 West Cummings Park, Suite 4950

Woburn, MA 01801

💻 781-376-0808

✉ inquiries@jalimedical.com

🌐 www.jalimedical.com

September 2009



US Distributor of Magstim Magnetic Stimulators

© 2009 Jali Medical, Inc.