

US Distributor of Magstim Magnetic Stimulators



# JALI MEDICAL

*introduces*

**magstim<sup>®</sup>**



## *Magstim 200<sup>2</sup> Magnetic Stimulator*

CREATING THE FUTURE THROUGH INNOVATION

Based on the original Magstim 200 design the new Magstim 200<sup>2</sup> monophasic magnetic nerve stimulator draws upon the years of knowledge and expertise gathered in both the clinical and research environments, resulting in an innovative design that will establish the new Magstim as the system preferred choice.

PLEASE NOTE: Magstim 200<sup>2</sup> is intended for the stimulation of peripheral nerves.

# Magstim 200<sup>2</sup> specifications

## Power Input:

110 - 120V, 50 - 60Hz

## Pulse Characteristics:

100µs approximate rise time, 1ms total duration

## Repetition Rate:

Once every 2 seconds at up to 0 - 49% power level

Once every 3 seconds at up to 50 - 79% power level

Once every 4 seconds at up to 80 - 100% power level

## Capacitor Life Expectancy:

2 x 10<sup>6</sup> discharges at 70% power level

8 x 10<sup>5</sup> discharges at 80% power level

4 x 10<sup>5</sup> discharges at 90% power level

2 x 10<sup>5</sup> discharges at 100% power level

## Remote Coil:

Allows complete control of the system, maneuvering the power level, arming and triggering the unit

## UI Controller:

Easy to read digital display of the power level, coil temperature and system status

## Compatibility:

Compatible to all existing Magstim coils

## Approx. Dimensions:

18.1" Wide x 6.3" High x 14.8" Deep

## Approx. Weight:

45 LB



JALI MEDICAL, INC. Inquiries: 888-959-JALI

✉ 500 West Cummings Park, Suite 4950

Woburn, MA 01801

☎ 781-376-0808

✉ inquiries@jalimedical.com

🌐 www.jalimedical.com



October 2008

US Distributor of Magstim Magnetic Stimulators

© 2008 Jali Medical, Inc.